

## Microlabequipment-ELECTRIC ICU BED manufacturer in india (MTHF 10EIB)



1. Advance Linear Actuator system 04 Pcs Motor, 01 Pc Control Box, 01 Pc Handset.
2. Four Pieces individual side rails, ABS Moulded.
3. ABS Moulded Head and Foot panel with locking and corner buffer.
4. Frame work made of Rectangular M.S. Tube.
5. 4 section top made up of Oval perforated CRCA MS Sheet.
6. Bed Mounted on 5" Dia Castor with central braking system.
7. S.S. Telescopic IV Rod with 4 Locations
8. Overall Approx Size: 2090L x 910W x 460-700H MM.
9. Finish: Pretreated & Epoxy Powder Coated
10. Positions obtained by hand Remote.

### **Description :**

1. Advance Linear Actuator system 04 Pcs Motor, 01 Pc Control Box, 01 Pc Handset.
2. Four Pieces individual side rails, ABS Moulded.
3. ABS Moulded Head and Foot panel with locking and corner buffer.
4. Frame work made of Rectangular M.S. Tube.
5. 4 section top made up of Oval perforated CRCA MS Sheet.
6. Bed Mounted on 5" Dia Castor with central braking system.
7. S.S. Telescopic IV Rod with 4 Locations
8. Overall Approx Size: 2090L x 910W x 460-700H MM.
9. Finish: Pretreated & Epoxy Powder Coated
10. Positions obtained by hand Remote.

### **Specifications :**

ELECTRIC ICU BED manufacturer in india

---

microlabequipment

[www.microlabequipment.com](http://www.microlabequipment.com) | Email: [info@microlabequipment.com](mailto:info@microlabequipment.com)