

Microlabequipment-Plain Beds with SS Panel manufacturer in india (MTHF 09PBAP)



- Frame work made M.S Tube.
- Top made of perforated MS Sheet.
- ABS head & foot bows of unequal height.
- Location for four I.V. Rod.
- Without Castor wheels..
- Overall Size: 2060 (L) x 900 (W) x 600 (H).

Description :

- Frame work made M.S Tube.
- Top made of perforated MS Sheet.
- ABS head & foot bows of unequal height.
- Location for four I.V. Rod.
- Without Castor wheels..
- Overall Size: 2060 (L) x 900 (W) x 600 (H).

Specifications :

Plain Beds with SS Panel manufacturer in india