

Microlabequipment-Smoke Density Apparatus manufacturer in india (MTEM-03)



Smoke Density Apparatus offered meets demands of measuring optical density of smoke upon certain thick specimen. Some of its features include:

- Widely used in assessing characteristics of vehicle interior and inflammable cables
- Integrated designs
- Auto open ventilation to control exhaust more efficiently
- Teflon coated inside for easy cleaning and corrosion resistant usage

Description :

Smoke Density Apparatus offered meets demands of measuring optical density of smoke upon certain thick specimen. Some of its features include:

- Widely used in assessing characteristics of vehicle interior and inflammable cables
- Integrated designs
- Auto open ventilation to control exhaust more efficiently
- Teflon coated inside for easy cleaning and corrosion resistant usage
- Digital display of smoke optical transmittance and density
- Pre-heated apparatus for quick start up and accurate data acquiring
- Safe blowout panel for safety
- Weighing sensor for enabling mass loss rate accurately

Specifications :

| | |
|--|--|
| | 25 KW/m ² or 50 KW/m ² |
| | Mild Steel |
| | 0 - 100% |
| | 300 x 300 x 790 mm |

230 V AC

microlabequipment

www.microlabequipment.com | Email: info@microlabequipment.com