

Microlabequipment-Wooden Lab Stools manufacturer in india (MTF 01WLS)



Wooden Lab Stools manufacturer in india

Description :

Wooden Lab Stools manufacturer in india

Specifications : Shape of the stool Top

Square

Style of Stools

Supported by set of legs

Loading bearing capacity

75 kg

Seat Support

Yes

Leg Support

All four Side

Footrest

According to requirement

Leg cross section members

Yes

Type of joints

Bridle joint

Mode of supply

Assembled ready to use

Top Made up of

Teak Wood

Seat support made up of

Teak Wood

Frame made up of

Teak Wood

Dimension of Table top	350 (L) mm x 350 (W) mm
Height of Stool ± 10 mm	600 mm
Thickness of Table top	25 mm (+/-1 mm)
Leg sizes	35 (L) mm x 35 (W) mm
Dimension of seat support	60 mm X 20 mm (thickness)
Leg support/cross section size	45 mm x 45 mm
Stool top finish	Lacquered polish
Frame finish	Lacquered polish

microlabequipment

www.microlabequipment.com | Email: info@microlabequipment.com