

Microlabequipment-2 Shelf Instrument trolley manufacturer in india (MTHF 04SIT)



- Stainless steel.
- Three shelves.
- Size: 760*460*810mm.
- With handle bar on short side.
- Noiseless casters with brake.

Description :

- Stainless steel.
- Three shelves.
- Size: 760*460*810mm.
- With handle bar on short side.
- Noiseless casters with brake.

Specifications :

Type of Instrument Trolley	2 Shelf
Shelves Position	one on
No of Caster	4
Availbility of Breaks	2 Cas
Availbility of Chrome Plated Pressed Sheet Frame	yes
Availbility of Raling above the Shelves	Yes
Space of Railing	50 mm
Diameter of Railing	8 mm

Thickness of Railing Sheet	6mm
Foldability	No
All Stainless Steel parts should be lightly polished in Dull Buff finish	yes
In order to avoid sharp edge the bent edge should be headed	yes
Chrome plated Frame work	yes
Quality of Paint on the Trolley must be in accordance with latest IS for paints	yes
Bowl Holder and a Stainless Steel Bowl and a Plastic Pail	Yes
Load carrying capacity	70 kg
Weight carrying capacity	30 kg
Mode of operation	Manu
Warranty	1 Yea
Material of trolley	Stainl
Thickness of Sheet of trolley in mm	6
Thickness of Sheet of Shelf in mm	1.6
No of Vertical Pole	4
Diameter of Vertical Pole in mm	30
Thickness of Sheet of Vertical Pole in mm	1.6
Length of Trolley in mm	760
Width of Trolley in mm	460
Height of Trolley in mm	810
Overall Weight of the trolley (Kg)	16
Material of Caster	Rubb
Diameter of Caster in mm	100
Ground Clearance in mm	150