

Microlabequipment-ADVANCED FULL FUNCTIONAL THREE- YEAR OLD CHILD NURSING MANIKIN manufacturer in india (MT-YM-215)



Features:

1. Bilateral Pupil observation: normal, mydriasis and myosis etc.
2. Operations of child care: full body bandage exercise, bathing, changing clothes, hair combing, eye washing and administering, otic drops and irrigation mouth cavity care, blood collection of finger and TB test.
3. Airway administering: realistic moth, nose, tongue, gum, pharynx, larynx, esophagus, epiglottis, trachea and tracheal ring, can practice tracheal intubation, sputum suction, oxygen inhaling

Description :

Features:

1. Bilateral Pupil observation: normal, mydriasis and myosis etc.
2. Operations of child care: full body bandage exercise, bathing, changing clothes, hair combing, eye washing and administering, otic drops and irrigation mouth cavity care, blood collection of finger and TB test.
3. Airway administering: realistic moth, nose, tongue, gum, pharynx, larynx, esophagus, epiglottis, trachea and tracheal ring, can practice tracheal intubation, sputum suction, oxygen inhaling
4. Oral and nasal feeding
5. Gastrolavage
6. Venipuncture, injection, transfusion (Arm). There is an obvious feeling when needle penetrates into vein and blood flash back indicates proper insertion.
7. Intramuscular injection in bilateral deltoid, bilateral vastus lateralis and bilateral buttock. Enema, placement of rectal suppository
8. Urethral catheterization
9. Ostomy nursing of ileum, colon and bladder.
10. Holistic care : sponge bath, replacing clothes, cold and heat therapy.
11. Limbs joints can bend, rotation, upper and lower movement.

12. Bathing and bandaging activity

13. Tibial bone marrow puncture: simulate the process of bone marrow puncture, obvious body surface symbol. Simulated bone marrow would outflow in the case of correct puncture.

Specifications :

ADVANCED FULL FUNCTIONAL THREE- YEAR OLD CHILD NURSING MANIKIN manufacturer in india

microlabequipment

www.microlabequipment.com | Email: info@microlabequipment.com